



iConnect[®] Enterprise Archive

V 11.4

IHE INTEGRATION STATEMENT

Merge Healthcare
900 Walnut Ridge Drive,
Hartland, WI 53029
USA

MERGE
An IBM[®] Company

877.44.MERGE • merge.com

 [@MergeHealthcare](https://twitter.com/MergeHealthcare)

 linkedin.com/company/merge-healthcare

 facebook.com/MergeHealthcare

© Copyright 2017 Merge Healthcare Incorporated, an IBM Company.

The content of this document is confidential information of Merge Healthcare Incorporated and its use and disclosure is subject to the terms of the agreement pursuant to which you obtained the software that accompanies the documentation.

Merge Healthcare® is a registered trademark of Merge Healthcare Inc.

The Merge Healthcare logo is a trademark of Merge Healthcare Inc.

All other names are trademarks or registered trademarks of their respective companies.

U.S. GOVERNMENT RESTRICTED RIGHTS:

This product is a “Commercial Item” offered with “Restricted Rights.” The Government's rights to use, modify, reproduce, release, perform, display or disclose this documentation are subject to the restrictions set forth in Federal Acquisition Regulation (“FAR”) 12.211 and 12.212 for civilian agencies and in DFARS 227.7202-3 for military agencies. Contractor is Merge Healthcare.

INTENDED USE:

Merge iConnect Enterprise Archive is intended for use as a vendor neutral archive for storage and communications of medical images and data.

The symbols glossary is provided electronically at <http://www.merge.com/Support/Resources.aspx>.

CAUTION: U.S. federal law restricts this device to sale by, or on the order of, a physician.

CANADIAN DEVICE IDENTIFIER: 40020



EMERGO EUROPE
Prinsessegracht 20
2514 AP, The Hague
The Netherlands

Australian Sponsor
IBM Australia Limited
Level 13, IBM Centre
601 Pacific Highway
St Leonards, NSW, 2065
Australia
ABN 79 000 024 733
Phone: 1800 117 425

Manufacturer's Address
Merge Healthcare Incorporated
900 Walnut Ridge Drive
Hartland, WI 53029

For application support or to report issues with user documentation, contact Customer Support:

- In North America: call toll-free 1-877-741-5369
- Outside of North America: +31.40.299.0773
- Email: MergeSupport@us.ibm.com

Part	Date	Version	Description
ECM-13050	July 2017	1.0	Initial version for 11.4

The latest revision of this document can be found in [Merge Communities](#).

Chapter 1 IHE Integration Statement

The following table outlines the IHE Integration Statement.

IHE Integration Statement		Date
		2015
Vendor	Product Name	Version
Merge Healthcare	iConnect® Enterprise Archive	11.1
This product is intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Patient Information Reconciliation	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Consistent Presentation of Images	Image Manager / Image Archive	None
Presentation of Grouped Procedures	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Evidence Documents	Image Manager / Image Archive	None
Key Image Note	Image Manager / Image Archive	None
Simple Image and Numeric Report	Report Repository	None
	Enterprise Report Repository	None
	External Report Repository Access	None

Access to Radiology Information	Image Manager / Image Archive	None
	Report Repository	None
	External Report Repository Access	None
Nuclear Medicine	Image Manager / Image Archive	None
Cardiac Catheterization Workflow	Image Manager / Image Archive	Intermittently Connected Modality Cardiac Cath
	Performed Procedure Step Manager	None
Echo Workflow	Image Manager / Image Archive	Echocardiography Intermittently Connected Modality
	Performed Procedure Step Manager	None
Evidence Documents - Cardiology	Image Manager / Image Archive	None
Consistent Time (CT)	Time Client	None
Mammography	Image Manager / Image Archive	None
Stress Testing Workflow	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Image Document Source	Set of DICOM Instances
	Document Repository	None
Cross-Enterprise Document Sharing (XDS.b)	Document Repository	None
	Integrated Document Source/Repository	Document Replace
Audit Trail and Node Authentication (ATNA)	Secure Application	None
Patient Identifier Cross-reference (PIX)	Patient Identifier Cross-reference Consumer	None
Patient Demographics Query (PDQ)	Patient Demographics Consumer	None

Imaging Object Change Management (IOCM)	Image Manager/Archive	None
	Change Requester	None
Internet address for vendor's IHE information:	www.merge.com/Support/Resources/IHE-Statements.aspx	
Links to Standards Conformance Statements for the Implementation		
HL7	www.merge.com/Support/Resources/HL7-Specs-and-Statements.aspx	
DICOM	www.merge.com/Support/Resources/DICOM-Conformance-Statements.aspx	
Links to general information on IHE		
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org