

| IHE Integration Statement |                 |         |        |
|---------------------------|-----------------|---------|--------|
| Vendor                    | Product Name    | Version | Date   |
| Merge Healthcare          | iConnect Access | 6.3     | 2/3/17 |

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

| Integration Profile                                    | Actor                                       | Integration Profile Option |
|--|---|----------------------------|
| Audit Trail and Node Authentication                    | Secure Application                          | None                       |
| Consistent Time  | Time Client                                 | None                       |
| Patient Demographics Query                             | Patient Demographics Consumer               | None                       |
| Patient Identifier Cross-referencing for MPI           | Patient Identifier Cross-reference Consumer | None                       |
| Cross-Community Access for Imaging                     | Imaging Document Consumer                   | None                       |
| Cross-Enterprise Document Sharing                      | Document Consumer                           | None                       |
| Cross-Enterprise Document Sharing for Imaging          | Imaging Document Consumer                   | None                       |
| Cross-Enterprise Document Sharing of Scanned Documents | Content Consumer                            | None                       |
| Cross-Enterprise Document Sharing of Scanned Documents | Content Consumer                            | View Option                |

Internet address for vendor's IHE information :  
<http://www.merge.com/Support/Resources/IHE-Statements.aspx>

Links to Standards Conformance Statements for the Implementation

|              |   |
|--------------|---|
| <b>HL7</b>   |   |
| <b>DICOM</b> | <a href="http://www.merge.com/MergeHealthcare/media/support/interop/iConnectAccessDCS6-3.pdf">http://www.merge.com/MergeHealthcare/media/support/interop/iConnectAccessDCS6-3.pdf</a> |

Link to general information on IHE : [www.ihe.net](http://www.ihe.net)