



IHE Integration Statement
Merge Eye Station™
v. 11.3

Merge Healthcare
900 Walnut Ridge Drive
Hartland, WI 53029
USA
877.44.MERGE

© 2012 Merge Healthcare. The information contained herein is confidential and is the sole property of Merge Healthcare. None of the information may be reproduced, copied, and/or re-distributed electronically, manually, or otherwise to any person without advance written permission from Merge Healthcare, except as may be expressly authorized by contracts or other agreements between Merge Healthcare and its customer. Merge Healthcare® and AMICAS are registered trademarks and VERICIS, RadSuite, Abraxas, and iConnect are also trademarks of Merge Healthcare. Other trademarks or service marks referred to herein are the property of their respective owners.

Merge Eye Station was formerly known as OIS WinStation.

Doc #	Revision	Date	Description
WSTN-9131	1	8/2/12	Initial release for 11.3

IHE Integration Statement		Date
		2012
Vendor	Product Name	Version
Merge Healthcare	Eye Station	11.3
This product is intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Eye Care Workflow	Acquisition Modality	PPS Exception Management
		Eye Care Image
Charge Posting	Acquisition Modality	None
<u>Internet address for vendor's IHE information:</u>		www.merge.com/Support/Resources/IHE-Statements.aspx
Links to Standards Conformance Statements for the Implementation		
HL7	www.merge.com/Support/Resources/HL7-Specs-and-Statements.aspx	
DICOM	www.merge.com/Support/Resources/DICOM-Conformance-Statements.aspx	
Links to general information on IHE		
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org