



IHE Integration Statement
iConnect[®]
Enterprise Archive
v. 9

Merge Healthcare
900 Walnut Ridge Drive
Hartland, WI 53029
USA
877.44.MERGE

© 2013 Merge Healthcare. The information contained herein is confidential and is the sole property of Merge Healthcare. None of the information may be reproduced, copied, and/or re-distributed electronically, manually, or otherwise to any person without advance written permission from Merge Healthcare, except as may be expressly authorized by contracts or other agreements between Merge Healthcare and its customer. Merge Healthcare® and AMICAS are registered trademarks and VERICIS, RadSuite, Abraxas, and iConnect are also trademarks of Merge Healthcare. Other trademarks or service marks referred to herein are the property of their respective owners.

Canadian Device 40020

**Australian Sponsor
(Regulatory Affairs Only)**
Emergo Asia-Pacific Pty Ltd.
Darling Park, Tower II, Level 20
Sydney, NSW 2000 Australia

Telephone numbers

Telephone: +61 (0) 2.9006.1662
Fax: +61 (0) 2.9006.1010



EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands
Tel: +31 (0)70 345 8570
Fax: +31 (0)70 346 7299

| Doc # | Revision | Date | Description |
|----------|----------|---------|------------------------------------|
| ECM-2070 | 2 | 2/20/13 | Addedd Canadian device identifier. |

| IHE Integration Statement | | Date |
|--|-----------------------------------|---|
| | | 2012 |
| Vendor | Product Name | Version |
| Merge Healthcare | iConnect® Enterprise Archive | 9 |
| This product is intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: | | |
| Integration Profiles Implemented | Actors Implemented | Options Implemented |
| Scheduled Workflow | Image Manager / Image Archive | None |
| | Performed Procedure Step Manager | None |
| Patient Information Reconciliation | Image Manager / Image Archive | None |
| | Performed Procedure Step Manager | None |
| Consistent Presentation of Images | Image Manager / Image Archive | None |
| Presentation of Grouped Procedures | Image Manager / Image Archive | None |
| | Performed Procedure Step Manager | None |
| Evidence Documents | Image Manager / Image Archive | None |
| Key Image Note | Image Manager / Image Archive | None |
| Simple Image and Numeric Report | Report Repository | None |
| | Enterprise Report Repository | None |
| | External Report Repository Access | None |
| Access to Radiology Information | Image Manager / Image Archive | Multiple Sources |
| | Report Repository | Multiple Sources |
| | External Report Repository Access | None |
| Nuclear Medicine | Image Manager / Image Archive | None |
| Cardiac Catheterization Workflow | Image Manager / Image Archive | Intermittently Connected Modality Cardiac Cath |
| | Performed Procedure Step Manager | None |
| Echo Workflow | Image Manager / Image Archive | Echocardiography Intermittently Connected Modality |
| | Performed Procedure Step Manager | None |
| Evidence Documents - Cardiology | Image Manager / Image Archive | None |
| Consistent Time (CT) | Time Client | None |
| Mammography | Image Manager / Image Archive | None |

| | | |
|--|--|--|
| Stress Testing Workflow | Image Manager / Image Archive | None |
| | Performed Procedure Step Manager | None |
| Cross-Enterprise Document Sharing for Imaging (XDS-I.b) | Imaging Document Source | Set of DICOM Instances |
| Audit Trail and Node Authentication (ATNA) | Secure Node | ATNA Encryption |
| Patient Identifier Cross-reference (PIX) | Patient Identifier Cross-reference Consumer | None |
| Patient Demographics Consumer (PDQ) | Patient Demographics Consumer | None |
| <u>Internet address for vendor's IHE information:</u> | | www.merge.com/Support/Resources/IHE-Statements.aspx |
| Links to Standards Conformance Statements for the Implementation | | |
| HL7 | www.merge.com/Support/Resources/HL7-Specs-and-Statements.aspx | |
| DICOM | www.merge.com/Support/Resources/DICOM-Conformance-Statements.aspx | |
| Links to general information on IHE | | |
| In North America: http://www.ihe.net | In Europe: http://www.ihe-europe.net | In Japan: http://www.ihe-j.org |